



# Product Specification Sheet

DSM Food Specialties B.V.

P.O. Box 1  
2600 MA Delft  
The Netherlands

[www.dsm.com](http://www.dsm.com)

## Delvotest® Fast Go Max

### Description

Delvotest® Fast Go Max (“Reader”), is a reader specifically designed to be used with DSM’s Delvotest® Fast range of tests for the detection of antibiotic residues in milk in cassette presentations, it gives you qualitative assessment by capturing and analyzing images of lateral flow immunochromatography.

### Technical data

<b>Model name</b>	Delvotest® Fast Go Max
<b>Operating System</b>	Embedded Linux
<b>Display</b>	10” Color LCD (Resolution 800 x 480)
<b>I/O Interface</b>	4 x USB 2.0 / 1 x RJ-45 LAN Connector
<b>Control</b>	Touch Screen
<b>Supported Accessories</b>	Barcode reader, Keyboard, Mouse, USB Memory (Optional)
<b>Input Voltage</b>	12 VDC
<b>Power Consumption</b>	12W
<b>Dimension (W x H x D)</b>	300 mm x 186 mm x 210 mm
<b>Weight</b>	1.8 kg
<b>Operating Temperature</b>	10° C ~ 40° C
<b>Storage Temperature</b>	-10° C ~ 50° C
<b>Features</b>	High sensitivity & repeatability 1D barcode support

Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty for which we assume legal responsibility, including without limitation any warranties as to the accuracy, currency or completeness of this information or of non-infringement of third party intellectual property rights. The content of this document is subject to change without further notice. This document is non-controlled and will not be automatically replaced when changed. Please contact us for the latest version of this document or for further information. Since the user’s product formulations, specific use applications and conditions of use are beyond our control, we make no warranty or representation regarding the results which may be obtained by the user. It shall be the responsibility of the user to determine the suitability of our products for the user’s specific purposes and the legal status for the user’s intended use of our products.

DSM Food Specialties B.V. | A. Fleminglaan 1 | 2613 AX Delft | The Netherlands | [info.food@dsm.com](mailto:info.food@dsm.com) | Trade Register Number 27235314



## Environmental

For indoor use only  
Main supply voltage fluctuations should not exceed  $\pm 10\%$   
Pollution degree 2

## Warranty:

Unless stipulated otherwise in the relevant sales agreement, DSM provides a guarantee on the Reader only, which guarantee is stipulated in the General terms and Conditions of Sale of DSM ("GTCS"). Supported Accessories are explicitly excluded from any guarantee given by DSM. The Reader is guaranteed to be free from defective materials and workmanship within the defined period, provided that misuse, alteration and/or abnormal use have not occurred by user. Do not open or repair defective readers but contact DSM representative immediately when a defect occurs. The Reader will be repaired under warranty for a period of 1 (one) year from the date of delivery of the device. Products NOT repaired by DSM (third party repair) are explicitly excluded from any guarantee or warranty given by DSM. In case of contradictions between these specifications and the GTCS, the GTCS shall prevail.

Although diligent care has been used to ensure that the information provided herein is accurate, nothing contained herein can be construed to imply any representation or warranty for which we assume legal responsibility, including without limitation any warranties as to the accuracy, currency or completeness of this information or of non-infringement of third party intellectual property rights. The content of this document is subject to change without further notice. This document is non-controlled and will not be automatically replaced when changed. Please contact us for the latest version of this document or for further information. Since the user's product formulations, specific use applications and conditions of use are beyond our control, we make no warranty or representation regarding the results which may be obtained by the user. It shall be the responsibility of the user to determine the suitability of our products for the user's specific purposes and the legal status for the user's intended use of our products.

DSM Food Specialties B.V. | A. Fleminglaan 1 | 2613 AX Delft | The Netherlands | [info.food@dsm.com](mailto:info.food@dsm.com) | Trade Register Number 27235314

Date of issue: February 18, 2019

© DSM Food Specialties B.V.

Page 2 of 2

HEALTH • NUTRITION • MATERIALS